

THD Basic Plus

Code 800075

20 sterile Kits

Each Kit contains:

- 1 THD device
- 1 THD Surgy Maxi
- 1 set of absorbable sutures

2 multi-use Doppler probes (10 shots each)

1 reusable needle holder



THD Basic

Code 80049C

20 sterile Kits

Each Kit contains:

- 1 THD device
- 1 THD Surgy Mini

2 multi-use Doppler probes (10 shots each)

1 reusable needle holder



Code 80049D

20 sterile Kits

Each Kit contains:

- 1 THD device
- 1 THD Surgy Maxi

2 multi-use Doppler probes (10 shots each)

1 reusable needle holder

Optional



THD Evolution

Code 800000



THD ShiningLight

Code 800051



THD S.p.A.
Via Industria, 1
42015 Correggio (RE) Italy

Customer Service
Ph. +39 0522 634343 - Fax +39 0522 693552
www.thdlab.com

360118 - 03/2011/0

THD S.p.A. reserves the right to change the products without prior notice.

THD Basic

The Smart Solution



THD Basic

The Smart Solution

DEARTERIALIZATION

Identification and selective ligation of the terminal branches of the superior rectal artery (6 vessels).

Standardised and Reliable

The guide into which the tip of the needle holder is inserted guarantees standardised and repeatable ligation and stitch depth control.



Doppler Guide



Precise ligation



NO excision



THD Basic



THD Basic

CAN BE USED TOGETHER WITH:
THD Surgy Maxi or THD Surgy Mini

MUCOPEXY

Repositioning of the cushions in their anatomical site by folding the mucohaemorrhoidal prolapse.

Functional Illumination

The design and the optical fibre system provides excellent illumination of the tissue being treated for both THD Basic and THD Surgy.

Optimal Surgical Field

THD Basic and THD Surgy have been designed to offer high functionality and to provide an ideal surgical field and perfect operating access. THD Basic device is optimized for dearterialization, while THD Surgy provides optimal exposure during pexy of the rectal mucosa.



THD Surgy Maxi



THD Surgy Mini