



An Intelligent Angle on Laparoscopic Fixation



Security in Hand

- Absorbable, spiral tacks provide secure mesh fixation
- Allows 360° fixation of the mesh from one side, reducing the need for placing additional trocars on the contralateral side
- An in-line handle provides a more ergonomic design that reduces stress and increases surgeon comfort





An Intelligent Angle on Laparoscopic Fixation

The iMESH Tacker™ incorporates a proprietary angling tip and in-line handle to enable a secure tack fixation angle through fewer trocars, leading to improved outcomes and optimal patient care.

Other tackers may require multiple trocar sites to securely tack the mesh, potentially leading to complications and longer operating time. The iMESH Tacker™ angulation allows fixation of the mesh 360° from one side, reducing the need for placing additional trocars on the contralateral side. The in-line handle provides a more ergonomic design that reduces stress and increases surgeon comfort.


The iMESH Tacker™ can be loaded in an angled position for intraoperative efficiency. The tacks are absorbable and longer than other available tacks. The additional length combined with the angulating tip provide the laparoscopic surgeon more consistent secure mesh fixation.

Together with the iMESH Placer™ and iMESH Stitcher™, the iMESH Tacker™ makes laparoscopic hernia repair more consistent and efficient than ever before.

The iMESH Placer™ and iMESH Stitcher™ are FDA cleared. The iMesh Tacker™ is awaiting FDA clearance.



An in-line handle provides a more ergonomic design that reduces stress and increases surgeon comfort.



An articulating tip enables reaching difficult positions.

Tacker Handle



The ergonomic, in-line handle design allows the surgeon to hold the instrument comfortably throughout the entire fixation process.

Tacker Tip



An articulating tip enables tacking from both the lateral and contralateral sides, and reaching difficult positions.

Tack



Absorbable, spiral tacks provide secure mesh fixation.



info@easy-lap.com www.easy-lap.com

Headquarters

30 Ha'Ella St.
POB 128 Kfar Truman, Israel 73150
+972 (77) 6201003

Easy-Lap Sales in the United States

190 Industrial Road, Suite 2
Wrentham, MA 02093, USA
(508) 384-8490

The iMESH Placer™, iMesh Tacker™ and iMESH Stitcher™ are patents pending.